Radiation Physics and Protection Series Units 1 – 12 DVD Format

Prepared by: Carolyn J. Frigmanski, M.A., B.S.R.T. [®] Founder, S.T.A.R.S.

All post tests must be returned for the designated number of Ohio Department of Health approved continuing education credits.



Specialized Topics in Areas of Radiologic Sciences

P.O. Box 2931 Toledo, Ohio 43606 419-471-1973

E-mail: info@xrayhomestudies.com

www.xrayhomestudies.com

Welcome to your **S.T.A.R.S.** DVD self-learning product that has been approved for **12 continuing education credit hours** by the Ohio Department of Health for radiologic license requirements. You may complete these units by using the audio-visual format from the Presentation Menu on the DVD or by printing the hard copy from the Resources Menu on your DVD. The total module is approved for **12 CE credits.** You must complete the post test in each unit with your name, address and date (social security number is optional) with a 75% or better score to receive an original signed certificate of completion with the respective CE credit amount.

Table of Contents

Unit 1	Radiologic Physics
Unit 2	Radiologic Physics continued with Basic X-ray Circuitry
Unit 3	Radiologic Quality/Technique – Factors affecting Density
Unit 4	Radiologic Quality/Technique – Factors affecting Contrast
Unit 5	Radiologic Quality/Technique – Factors affecting Detail & Distortion
Unit 6	Radiologic Quality/Technique – Accessories for improving the Diagnostic Image
Unit 7	Darkroom Processing and Film Handling
	Digital Imaging Principles
Unit 8	Radiobiology – Cellular and Molecular Biology
Unit 9	Radiobiology – Acute Effects of Radiation Exposure
Unit 10	Radiobiology - Chronic Effects of Radiation Exposure
Unit 11	Radiation Protection for the Patient
Unit 12	Radiation Protection for the Operator

Please submit the post tests to the **S.T.A.R.S.** office **BEFORE** your biennium ends to receive CE credit. Certificates of completion are prepared with the date your post tests are received and evaluated by **S.T.A.R.S.** officials.

Thank you very much for selecting this educational product to meet your professional needs.

Carolyn J. Frigmanski, M.A.,B.S.R.T.® Founder



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo, Ohio 43606 419-471-1973

E-mail: info@xrayhomestudies.com

www.xrayhomestudies.com

Dear GXMO/LSO/LSR participant,

Welcome to your selected S.T.A.R.S. developed continuing education home study!

In the spring of 2013, The Ohio Department of Health (ODH) approved ALL of them for GXMO license renewal. Please check with your state licensing agency if you are not a licensed GXMO in Ohio to be sure your state accepts these ce credits for your state license renewal requirement.

This product consists of a text on a desired topic and multiple question, short answer post test(s) pages. The number of Ohio Department of Health approved continuing education (ce) credits is listed on our order form. This home study product was originally developed for radiographers across the country who were registered with the American Registry of Radiologic Technologists (ARRT) and approved by the American Society of Radiologic Technologists (ASRT). Please disregard any reference to the ASRT/ARRT within this product, if any.

You must complete the reading and questions with a **75% or higher score** on the post test(s) to get your approved CE credit!

Please return all the post test pages to the **S.T.A.R.S.** address listed on our letterhead **BEFORE** your license expires. We will forward your certificate of completion on the same day your post test(s) were evaluated (except for holidays and Sundays). If you did NOT receive a 75%, we will send the pages back with the questions needing a new answer selection. Send them back for a **FREE** reevaluation. No refund will be provided for unsatisfactory personal performance on any ce product.

Plan the return of your post test(s) pages in a timely manner. I cannot accept emailed or faxed copies since I need to retain my **ORIGINAL** records for the ODH for 3 years in case you may be audited.

Remember to get your envelope weighed at the post office whenever submitting more than 4 pages. The post office will return it to you if you have insufficient postage, thereby delaying it for my evaluation and your certificate of completion.

Feel free to contact me by email: info@xrayhomestudies.com or telephone: 419 471-1973 if you have any questions. Please share with others in the future.

Thank you for selecting **S.T.A.R.S.** to meet your continuing education needs!

Sincerely,



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo. Ohio 43606 419-471-1973

P.O. Box 2931 Toledo, Ohio 43606 E-mail: info@xrayhomestudies.com

www.xrayhomestudies.com

How to renew your GXMO license in Ohio:

The Ohio Dept. of Health (ODH) requires a minimum of 12 continuing education credits (ceus) to be completed every 2 years (your biennium) BEFORE your license expires. You may do more than 12 ceus, but not less than 12 ceus, if you so choose. Ceus in excess of 12 cannot be carried over to the next biennium.

You will receive a hard copy renewal notice by mail from the ODH 60 days **BEFORE** your license expires. It is your responsibility to amend your personal information to the ODH whenever you change your name, address or place of employment as soon as possible by using the ODH website or contacting the ODH by telephone at 614 752-4319 for assistance. Failure to receive an ODH notice is not an acceptable reason for failing to renew on time. You can add completion of clinical modules to your GXMO license on the ODH web site.

Your ODH notice informs you that you may renew online or request a hard copy form from them. You must have your S.T.A.R.S. certificate(s) of completion immediately available when you renew since your course title(s), number of ceus, and ODH accreditation number(s) and date(s) of completion are printed on it.

You can renew immediately when you receive your notice or you have 30 days to complete the renewal process and payment to the ODH after your license expires. Online renewal requires your credit card for payment. If you chose hard copy renewal, you may submit a check or money order.

You and/or your employer can view and/or print your renewed license on line upon completion of the process. Problems that exist with renewal should be addressed to the ODH by calling for assistance.

S.T.A.R.S. personnel are **NOT** responsible for your renewal. Please direct any questions or needed assistance with renewal to the ODH personnel.

GXMOs must notify the ODH in writing within 30 days of any changes in the physician providing direct supervision. If your scope of practice changes (e.g. chiropractic to podiatric), a competency form must be completed and submitted to the ODH.

You may also want to check the ODH web site periodically for changes that may have occurred during your biennium and to share this information with your co-workers and/or administrative staff members.

The ODH website is: http://www.odh.ohio.gov/odhPrograms/rp/rlic/ristatus.aspx

Email is: BRadiation@odh.ohio.gov

Thank you very much.



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo, Ohio 43606 419-471-1973

E-mail: info@xrayhomestudies.com

www.xrayhomestudies.com

Instructions for Mailing your Continuing Education Post Tests

Complete ALL hard copy unit post tests for the products you purchased in legible printing BEFORE your license expires. Mail is processed the same day it is received.

You may want to copy them BEFORE you mail them to the S.T.A.R.S. office to minimize mail delivery complications. They will NOT be returned to you unless you get a 75% or less. If you do NOT get a 75% or better after evaluation, the post tests will be sent back to you with the questions needing a new answer selection. After completing the questions, send them ALL back to the S.T.A.R.S. office for re-evaluation.

Be sure to use the CORRECT postage by having it weighed at the post office if it consists of more than 5 pages. Envelopes with INSUFFICIENT POSTAGE will be sent back to the participant and delay your post test evaluation and certificate creation.

I do NOT accept faxes since faxes fade over time and I need to keep my records for 3 years in case you would get audited by the Ohio Dept. of Health for some reason.

I do NOT accept scanned pages because I do NOT want you putting your private, personal information on the internet. I do NOT open attachments due to the threat of virus contamination that may jeopardize my web site and computerized databases.

Do NOT send your study media i.e. CD, DVD, booklets and/or books back to me.

United States Postal Service (USPS):

If you are using USPS for priority or express mailing, please keep your receipt with the tracking number in case of a problem with the delivery. Please mark the section for NO SIGNATURE REQUIRED for express mail and send it to my home address: Carolyn J. Frigmanski, MA, BSRT (R) 3134 Aldringham Road Toledo, Ohio 43606. The USPS does NOT deliver to my P.O. Box address. Please call to let me know I should be expecting it at 419 471-1973.

FedEx or United Parcel Service UPS:

If you are using these delivery services, please keep your receipts with the tracking number in case of a problem with the delivery. Please mark the section for NO SIGNATURE REQUIRED for express mail and send it to my home address: Carolyn J. Frigmanski, MA, BSRT(R) 3134 Aldringham Road Toledo, Ohio 43606. Please call to let me know I should be expecting it at 419 471-1973.

Thank you very much.



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo, Ohio 43606 419-471-1973

E-mail: info@xrayhomestudies.com

www.xrayhomestudies.com

Dear S.T.A.R.S. participant,

Welcome to my DVD on Radiation Protection that non-licensed applicants may use for preparation for the radiologic license examination or for continuing education by licensed gxmos! There are over 800+ slides in its creation that correlate to the didactic curriculum developed by the Ohio Department of Health for license applicants and my CE product. I wanted you to know that you can play this DVD on your computer to view the entire presentation in addition to printing selected pages or all of it! You can also play this DVD on your TVs too so that others can view it collectively. You can stop, start and pause whenever necessary. It must be used with the hard copy S.T.A.R.S. Study Guide materials for license exam success. No certificate will be awarded for watching its contents. You must successful complete and return all the post tests to the S.T.A.R.S. office for evaluation and creation of an official certificate to be sent to you.

The **Presentation Menu** contains the audio-visual content of each individual topic for exam and CE content that has been approved by the Ohio Department of Health. Please disregard the Unit 00 due to ODH rule changes in 2012. The **Resources Menu** allows you hard copy printing of selected pages or all pages in any topic or all of the entire documents. Please follow the directions on the DVD label for printing capability.

As an educator, I realize that you may become overwhelmed by the course content. The primary motivation for this DVD is providing you an opportunity to review the contents at your own leisure and in the convenience of your home or workplace. Repetition is the best form of learning so I hope you view it frequently as you prepare for the license examination or when submitting your CE post tests for evaluation.

This DVD will be edited as changes occur within the state law or as personal needs change for general x-ray machine operators in their endeavors to successfully pass the state exam on the first attempt or CE post test. Some of the charts, graphs and/or tables were retyped for better readability when printed. There is a 5 second pause between frames to allow note taking and comprehension and a 10 second pause for frames with greater content. You may perceive a slight difference in my voice. I had to switch recording studio developers during the production process so the sound fluctuates minimally.

If you encounter any difficulty with your DVD, please call Stuart, my web designer, who is listed on the back of the DVD case. He will assist you with your problem.

Please feel free to contact me if you have any suggestions for improvements. It is always nice to get feedback from you as the learner rather than me as the developer.

Best wishes for your success in your career development.



Specialized Topics in Areas of Radiologic Sciences

P.O. Box 2931 Toledo, Ohio 43606 E-mail: info@xrayhomestudies.com

www.xrayhomestudies.com

The times for the various units are shown below.

419-471-1973

Unit 0: 0:17:33

Unit 1: 1:07:32

Unit 2: 0:51:58

Unit 3: 0:33:57

Unit 4: 0:27:03

Unit 5: 0:22:03

Unit 6: 0:28:39

Unit 7: 1:14:59

Unit 8: 1:04:56

Unit 9: 0:18:33

Unit 10: 0:31:17

Unit 11: 0:29:25

Unit 12: 0:41:23

Unit Dt: 1:04:14

Stuart

Stuart Feldman

The Solar Stop

419-THE-STOP

Each unit of exam material has been timed for your convenience in allocating enough time to study without interruption of content. Don't forget to disregard Unit 0 since it was orientation and no longer applies to the 2012 ODH rule changes.

The numbers indicate the hours, minutes and seconds of dictation when the unit is not interrupted with stops, pauses and/or restarts. For example, Unit 1 is 1 hour: 7 minutes: 32 seconds in length.

Stuart is my web site and DVD designer. Feel free to contact him if you have any problem using your DVD at 419 843-7867.



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo, Ohio 43606 419-471-1973

E-mail: info@xrayhomestudies.com

www.xrayhomestudies.com

- ❖ Please do **NOT** return the entire product! It's yours to keep except for the post test pages! I hope it will be a useful resource for you and others in your practice.
- ❖ You only need to fill out the information re: your name as you want it to appear on your certificate, your desired number i.e. social security number, last 4 digits of s.s. number or your license number; your address, etc. if the product you selected has multiple pages. You can initial the other pages and secure the group of pages with a staple so they don't become separated.
- Please use the correct postage when you return your post tests so delays are not incurred from insufficient postage.
- ❖ You can do the selected product only **ONCE**. Make other choices in the future.
- ❖ Others can share your media by ordering additional post tests at a reduced price. Group discounts are available too. Go to the **GXMO Order Form** tab for details.
- ❖ Don't forget to use the **GXMO Order Form** tab to select your next self/distance learning product for your license renewal in the future.



Unit/s: 1

Category: Radiologic Physics

Completing these worksheets is an example of	energy!
2. The positive particle found in the nucleus is called a	
3. A indicates the number of electrons in the outermo	ost shell.
4. An(A) is an atom that has gained or lost an elec	ctron.
5. Identify 1 important contribution from the process of ionization below	:
6. The greatest source of radiation exposure on average to the U.S. popula	ation is
7 particles are composed of the He nucleus.	
8. "Hard" x-rays possess frequency &	wavelengths.
9. Ionizing radiation is located above on the ele	ectromagnetic spectrum.
10. "Can Gramma really use Val's initial R" represents the order of the El really represents the EMR called	MR spectrum. The word
11 rays originate by disturbing electrons and rays or	originate from the nucleus.
12. X-rays can produce and	changes in the human body.
13. Identify 2 properties of x-rays:	
1.	
2.	
14. The interaction that occurs when incident x-rays ionize atoms and chenergy is calledeffect.	ange direction with a loss of
15. The radiation that remains after attenuation to reach the image receptor	or is called
16. As x-ray energy increases, we will have photo	ns transmitted through tissue
17. The traditional unit of measurement for patient doses is the	
18 The metric unit of radiation used as an occupational dose equivalent	is called



Fill in the following information so that an official from S.T.A.R.S. when all			
Name			
Address			
Social Security Number (optional)	-	_	Date



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo, Ohio 43606 419-471-1973

P.O. Box 2931 Toledo, Ohio 43606 E-mail: info@xrayhomestudies.com www.xrayhomestudies.com

TT 1// A	0-4	Dadislania	Dlarentas
Unit/s: 2	Category:	Radiologic	Physics

·	
1. The process of electrification by rubbing object	s together is called
2. Materials that inhibit the flow of electrons in a c	conductor are called
3. Identify 1 law of electrostatics:	
4. The electrostatic force is	proportional to the product of the charges.
5. Electromotive force is another term for the	in our circuits.
6. An example of a circuit is whe	en 1 Christmas tree bulb goes out, they all go out!
7 current involves	s the oscillation of electrons in opposite directions.
8. General x-ray tube circuits utilize	mA and kV.
9. The is the negative terminal	contained in the glass envelope of the x-ray tube.
10. The target interaction involving the deceleration	n of electrons from the cathode is called
	radiation.
11. Identify 1 method of prolonging x-ray tube life	e:
12. Thefocal spot is the surface of	on the anode that encounters the electron stream.
13. An increase in kVp will	the penetrating ability of the x-ray beam.
14. The primary function of filtration is to eliminate	te wavelengths from the x-ray beam.
15. Total filtration in diagnostic x-ray tubes is equ	nivalent to mm Al.
16. The protective housing reduces leakage radiati	on to less than mR/hr. at 1 meter.
17. AEC can be constructed with a photomultiplie	er tube or an chambe



Specialized Topics in Areas of Radiologic Sciences 419-471-1973

P.O. Box 2931 Toledo, Ohio 43606 E-mail: info@xrayhomestudies.com www.xrayhomestudies.com

Fill in the following information so that you can receive a signed certificate of c	ompletion from
an official from S.T.A.R.S. when all the sections of the unit are completed.	

Address		



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo, Ohio 43606 419-471-1973

E-mail: info@xrayhomestudies.com www.xrayhomestudies.com

Unit/s: 3 Category: Radiographic Quality/Technique

1.	Exactness of anatomic structures within a useful density range describes the property called
2.	Unwanted/unnecessary optical densities on a radiographic image is called
3.	The prime factor of exposure that determines the number of x-rays is
4.	The term is the degree of overall blackness on the image.
5.	Higher mA selection will provide density.
6.	Faster film & screen speeds will density.
7.	Shorter S.I.D/F.F.D. will provide density.
8.	Increased filtration will x-ray quantity.
9.	The law that allows us to change mA and time to achieve the SAME density when adjustments an needed due to the patient or equipment limitations is called the
	law.
10	Calculate the mAs: 200 mA @ 3/10 or .3 seconds =mAs.
11	The law that we can utilize to reduce our occupational dose is called the
	law.
12	The law that allows us to manipulate techniques to maintain the same density when distance
	changes is called thelaw or formula.
13	Materials that do NOT allow x-rays to pass through are called
14	The most dense tissue in the human body is
	Osteoporosis is considered/classified as a pathology.
	The body habitus of emaciated/extremely thin patients is described as
17	The sensitometric/H & D/characteristic curve that contains the blackest density is
	The useful diagnostic range for image interpretation is to .



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo, Ohio 43606 419-471-1973

E-mail: info@xrayhomestudies.com www.xrayhomestudies.com

	nformation so that you can a A.R.S. when all the sections		npletion from
Name			

Social Security Number (optional)	_	_	Date	
Address				



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo, Ohio 43606 419-471-1973

E-mail: info@xrayhomestudies.com www.xrayhomestudies.com

Unit/s: 4	Category: Radiographic	c Quality/Techniqu	e
1. Radiograph	nic images possessing many shades of gray are de	escribed as	scale.
2. Radiograp	hic contrast is created/composed by both	and	contrast.
3. High or sho	ort scale contrast is produced by using		kVp.
4. Processing contrast.	g solutions that are too hot or incorrect will produ	ice images with	scale
5. Patients w	ith high atomic numbers will produce images wit	th	scale contrast.
6. Tissues wit	th greater mass density produce	scale of cont	rast.
7. The most in	mportant exposure factor that controls contrast &	e penetration of the p	oart is
8. Increasing	kilovoltage will result in	fog due to	o increased scatter.
9. The film p	roperty that allows minor errors with technique is	s called	•
10. The optim	al kVp range for small extremities is & A	AP/ lateral vertebral	column is
11	filtration represents the attenuation/abso	orption by the windo	w of the x-ray tube.
12. A	% increase in kVp will	the density on the f	inished radiograph.
13. Name 1 de	evice/method to control scatter creation:		

14. Scatter and patient dose will ______ when a larger x-ray beam field is used.

15. When exposure time increases, quality _____ and quantity _____.



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 419-471-1973

Toledo, Ohio 43606

E-mail: info@xrayhomestudies.com www.xrayhomestudies.com

Fill in the following information so that you can receive a signed certificate of completion from an official from S.T.A.R.S. when all the sections of the unit are completed.

Address		



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo, Ohio 43606 419-471-1973

P.O. Box 2931 Toledo, Ohio 43606 E-mail: info@xrayhomestudies.com www.xrayhomestudies.com

Unit/s: 5 Category: Radiographic Quality/Technique

1.	1. The ability to image small objects that have hig	h subject contrast is called
2.	2. The term refers t	
3.	3. The focal spot is project	cted to the patient and affects image sharpness.
4.	4. A focal spot provides	the best radiographic detail.
5.	5 SID/FFD provides the best	radiographic detail, but also increases patient dose
6.	6OID/OFD provides the best	st radiographic detail.
7.	7. The biggest detriment to detail is	motion.
8.	8. Identify 2 immobilization techniques commonl	y used in radiography: 1.
	and 2	•
9.	9. Identify1 communication skill to use with our p	patients:
10	10. Technical errors will create	distortion.
11	11. Distortion can be created by the misalignment	of either the (3 items): 1.
	2 and/or _	•
12	12. The process of increasing the size of an object	to improve diagnostic quality is called
13	13. The 4 factors of radiographic quality are	discussed in Unit 3;
	discussed in Unit 4 and	d and
*	discussed in Unit 5	



Specialized Topics in Areas of Radiologic Sciences 419-471-1973

P.O. Box 2931 Toledo, Ohio 43606 E-mail: info@xrayhomestudies.com www.xrayhomestudies.com

Fill in the following information so that you can receive a signed certificate	of completion from
an official from S.T.A.R.S. when all the sections of the unit are completed.	

Name				
Address				
Social Secur	ity Number (optional)	_	Date	



Specialized Topics in Areas of Radiologic Sciences Toledo, Ohio 43606 P.O. Box 2931 419-471-1973

E-mail: info@xrayhomestudies.com www.xrayhomestudies.com

Unit/s: 6 Category: Conclusion of Radiographic Quality/Technique

1.	The primary function of a grid is to
2.	A grid has canted lead strips.
3.	is the ratio of transmitted primary to secondary radiation through the grid
4.	The height of the grid strips divided by the width of the interspace defines grid
5.	Identify 1 problem using a grid that will result in poor image quality:
6.	Improper use of grids will produce radiographic images demonstrating
7.	Identify 1 factor to consider when using a grid correctly:
8.	A 10/12:1 grid ratio requires x more exposure than non-grid.
9.	When converting from non-grid to grid exposure correction, you must
1(O. Identify 1 application/reason to use a grid:
11	1. Patient doses always go up when ratio grids are used.
12	2. 2 techniques used to improve image quality when filming body parts of varying thicknesses are
	1
	2
13	3. Thinner body parts should be placed under the terminal of the x-ray tube.



Fill in the following information so that you can receive a signed certificate of completion from an official from S.T.A.R.S. when all the sections of the unit are completed.

Name				· · · · · · · · · · · · · · · · · · ·	
Address					
Social Secur	rity Numbe	r (optional)		Date	



Unit/s: 7 Category: Darkroom Processing & Film Handling	
1. The photoreceptive layer of a radiographic film is also called its	
2. During exposure, the silver portion of the emulsion crystal becomes	_charged.
3. The property that describes the film's ability to respond to light and/or radiation is called	
or	
4. The conversion of the latent image into a manifest image occurs during	
in the processing cycle.	
5. The straight line portion of the sensitometric curve is also called	
6. Identify 2 safety precautions when handling processing chemicals:	
1.	_
2.	
7. Identify 1 function of the fixing solution:	
8. Describe 1 method to preserve the developer strength:	
9. Any undesirable density on a radiographic image is called	
10. The purpose of replenishment in the automatic processor is to	
<u> </u>	
11. Dirt or chemical stains on rollers may produce an artifact described as	lines
12. Proper storage of radiographic film is % humidity.	
13. The emission of light in intensifying screens is the principle of	
14. The random interaction of x-rays with screens is	
15. Intensifying screens with large crystals in thick active layers will pa	
16. Identify 1 precaution to utilize when cleaning and caring for intensifying screens:	



Specialized Topics in Areas of Radiologic Sciences 419-471-1973

P.O. Box 2931 Toledo, Ohio 43606 E-mail: info@xrayhomestudies.com www.xrayhomestudies.com

Social Security Number (optional) _____-

Fill in the following information so that you can receive a signed certificate of completion from an official from S.T.A.R.S. when all the sections of the unit are completed.				
Name				
Address				
	· · · · · · · · · · · · · · · · · · ·			

Date



Unit/s: Digital Imaging	Category: Digital In	naging Principles
1. A	signal consists of a string of	pinary numbers called
2. The smallest component o	f the matrix is called the	
3. Increasing the window lev	vel	the image brightness.
4. Window level controls	& wind	low width controls
		nanufacturers to communicate is called
		X
7. The 3 steps to produce a co	ligital image include: 1	-
2	and 3	
		•
9	conversion DR utili	zes a-Selenium as an x-ray photoconductor
10. The device made of photo	sensitive silicon chips is calle	ed a
11. Identify 1 difference betw	een CR and Conventional film	n:
12. Identify 1 difference betw	een CR and DR:	
		called a
14. The "" number	reflects the center of the usal	le portion of the histogram.
		to a radiographic image is called
16. The postprocessing techni	ique that removes superimpos	



Specialized Topics in Areas of Radiologic Sciences 419-471-1973

P.O. Box 2931 Toledo, Ohio 43606 E-mail: info@xrayhomestudies.com www.xrayhomestudies.com

Fill in the	following infor	mation so that yo	ou can receive a	a signed cer	tificate of c	ompletion from	n
an official	from S.T.A.R.	S. when all the se	ections of the u	nit are com	pleted.		

Name			
Address _			
Social Sec	urity Number (optional	 Date	



ry:	Radiobiology
	ry:

1. A effect occurs r	nany months or years after the initi	ial radiation exposure.
2. The most prevalent atom in the human bod	ly is	•
3. The most radiosensitive macromolecule su	bject to radiation damage is	
4	act as enzymes, hormones &	antibodies.
5. The double helix structure describes	located in the	of the cell.
6. Genetic change is the result of a radiation in	nteraction called	
7. Aeffect occ	curs when the ionizing radiation hi	ts the target molecule.
8. Breakage of the thread or backbone of a m	nacromolecule is called	
9. Identify 1 observable effect of radiation in		
10	are organelles that produce energy	
11. Somatic cells divide by	; sex cells undergo	
12 phase is the rep	plication of parent DNA into 2 daug	ghter DNA molecules.
13. The most radiosensitive sex cell stage for	males is	& females
14. A classic example of intracellular repair/r	repopulation that most of us have ex	
15. Our most radiosensitive cell type is		: most radioresistant is
	<u>:</u>	
16 rems can cause deat	th to human beings.	
17. The law that describes the principles of co	ell radiosensitivity was developed b	by 2 Frenchmen named
	and	•
18 effects occur	r when given enough dose.	

19. Radiation protection guidelines are based on a	•
dose-response relationship.	
20. Increasing LET and oxygen concentration	biologic damage.
21. The principal radiation interaction in the human body i	involves the radiolysis of
22. A is a highly reactive & unstabl	le molecule that creates damaging products
23. The most radiosensitive period in our human lives occu	
24. A fraction of exposure to the pregnant woman will exp	
creation of	radiation.
25. Identify 2 methods of protecting the pregnant patient fine.1.2.	
Fill in the following information so that you can receive a sofficial from S.T.A.R.S. when all the sections are complete	signed certificate of completion from an ed.
Name	
Address	
Social Security Number (optional)	Date



Unit/s: 9 Category: Radiobiology

1.	Cell counts can be depressed with rads of x-ray exposure.
2.	radiation syndrome occurs when ionizing radiation occurs over a short
	time.
3.	Identify 2 early effects of radiation exposure:
	1 and 2
4.	The stage involves hemorrhage, dehydration and infection.
5.	The gastrointestinal syndrome occurs when radiation exposure exceeds rads.
6.	The 50 in LD 50/30 indicates the
Ur	nit 10 Category: Radiobiology
1.	X-rays were discovered by and natural radioactivity was
	discovered by and
3.	The first fatality from man-made radiation in the U. S. was
3.	Identify 1 radiation injury sustained by our early investigators:
4.	What historic occupational population experienced cataract formation?
5.	Identify 1 effect of radiation on skin:
6.	A significant number of children treated with radiation therapy for thymus enlargement developed
	cancers.
7.	What historic occupational population experienced bone cancer?
8.	What historic population experienced an increase in breast cancer?
9.	What historic occupational population experienced lung cancer?
10.	. Permanent sterility in both sexes can result with rads or Gys.

11	risk describes the late	radiation effects in large populations when the
	precise dose is unkno	wn.
12. Life span shortening	for radiation workers is	days of expected life for every 1 rad based
	on mouse research.	
13. A controversial study	on life-span shortening was	experienced by
		* .
14. The annual risk of dy	ing this year from 100 mrad	s of radiation is
14. The annual risk of dy	ing this year from 100 mrad	s of radiation is
4. The annual risk of dy Space for notes: DVD Fill in the following inform	ing this year from 100 mrad	s of radiation is & McKinley pictures! Sorry about error. ve a signed certificate of completion from an
4. The annual risk of dy Space for notes: DVD	ing this year from 100 mrad has misdictation on Edisor mation so that you can rece	& McKinley pictures! Sorry about error. ve a signed certificate of completion from an eted.
Fill in the following informatical from S.T.A.R.S.	ing this year from 100 mrad has misdictation on Edisor mation so that you can rece when all sections are compl	& McKinley pictures! Sorry about error. ve a signed certificate of completion from an eted.



Unit/s: 11 **Category: Radiation Protection**

1. The measurement & testing of radiographic equipment is called
2. The entire rack assembly should be cleaned
3. Developer, wash water and replenishment rates should be checked
4. Collimator field size and centering should be checked
5. Identify 1 method/principle we can use to reduce unnecessary patient doses below:
6. Identify 1 way to minimize positioning errors:
7. Identify 1 factor to consider for time selection:
8. Identify 1 consideration when using a technique chart to establish accurate exposure calculation.
9. The technique chart utilizing specific kVp for a specific exam is called
10. Identify 1 factor that is selected when constructing an automated exposure chart.
11. Correct chamber selection and patient centering is critical when using equipment
12. Gonad shielding should be used if the gonads lie within orcms. to the primary x-ray field
13. Which gonadal shielding policy is preferred and explain why:
14. Identify 1 method/principle to provide patient protection:
15. Identify the 2 most common non-ionizing imaging modality options for pregnant patients.
1 and 2
16. The federal law about privacy practices, legal duties & patient rights is called
17. Identify 2 necessary permanent markers on our radiographic images:



Specialized Topics in Areas of Radiologic Sciences 419-471-1973

P.O. Box 2931 Toledo, Ohio 43606 E-mail: info@xrayhomestudies.com www.xrayhomestudies.com

Fill in the following information so that you can receive a signed certificate of completion from
an official from S.T.A.R.S. when all the sections of the unit are completed.

Name				
Address				
	-	- \		
Social Security Number (optional)	_	_	Date	



Unit/s: 12 **Category: Radiation Protection**

1.	The 3 cardinal principles of radiation protection are,	and and
2.	ALARA stands for	
	Identify 1 concept of ALARA:	
	The monthly dose limit during pregnancy ismSv and	
	The cumulative whole-body dose limit for a 26 year old radiation worker is	
6.	Identify 2 personnel protective devices: 1	and
	2.	
7.	Anyone likely to exceed% of the annual limit must be monitored	
8.	The radiation monitoring device that utilizes a lithium fluoride chip is called a	
	dosimeter.	
9.	Identify 1 exposure dose that is provided on a personnel monitoring report belo	w:
10.	No one is allowed to wear the monitor.	
11.	What do you do if your personnel monitor is lost?	
12.	Identify 1 equipment design that protects us from radiation exposure:	
	A protective barrier is constructed of "Pb or _	
	The best method to immobilize patients is to use	÷
	What is the first step for pregnant radiation workers to do?	
16.	Identify 3 important aspects of radiation protection other than the 3 cardinal print ALARA: 1.	



Fill in the following information so that you can receive a signed certificate of completion from an official from S.T.A.R.S. when all the sections of the unit are completed.							
Name							
Address							
				. '			

Social Security Number (optional) _____-___-Date