

Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo, Ohio 43606

Phone: 419-471-1973

Website: www.xrayhomestudies.com

Unit 36

Radiographic Imaging & Exposure

by

Terri L. Fauber

Third Edition

Prepared by: Carolyn J. Frigmanski, M.A., B.S.R.T. ® Founder, S.T.A.R.S.



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo, Ohio 43606 419-471-1973

E-mail: info@xrayhomestudies.com

www.xrayhomestudies.com

Dear GXMO/LSO/LSR participant,

Welcome to your selected S.T.A.R.S. developed continuing education home study!

In the spring of 2013, The Ohio Department of Health (ODH) approved ALL of them for GXMO license renewal. Please check with your state licensing agency if you are not a licensed GXMO in Ohio to be sure your state accepts these ce credits for your state license renewal requirement.

This product consists of a text on a desired topic and multiple question, short answer post test(s) pages. The number of Ohio Department of Health approved continuing education (ce) credits is listed on our order form. This home study product was originally developed for radiographers across the country who were registered with the American Registry of Radiologic Technologists (ARRT) and approved by the American Society of Radiologic Technologists (ASRT). Please disregard any reference to the ASRT/ARRT within this product, if any.

You must complete the reading and questions with a **75% or higher score** on the post test(s) to get your approved CE credit!

Please return all the post test pages to the **S.T.A.R.S.** address listed on our letterhead **BEFORE** your license expires. We will forward your certificate of completion on the same day your post test(s) were evaluated (except for holidays and Sundays). If you did NOT receive a 75%, we will send the pages back with the questions needing a new answer selection. Send them back for a **FREE** reevaluation. No refund will be provided for unsatisfactory personal performance on any ce product.

Plan the return of your post test(s) pages in a timely manner. I cannot accept emailed or faxed copies since I need to retain my **ORIGINAL** records for the ODH for 3 years in case you may be audited.

Remember to get your envelope weighed at the post office whenever submitting more than 4 pages. The post office will return it to you if you have insufficient postage, thereby delaying it for my evaluation and your certificate of completion.

Feel free to contact me by email: info@xrayhomestudies.com or telephone: 419 471-1973 if you have any questions. Please share with others in the future.

Thank you for selecting **S.T.A.R.S.** to meet your continuing education needs! Sincerely,



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo, Ohio 43606 419-471-1973

E-mail: info@xrayhomestudies.com

www.xrayhomestudies.com

Instructions for Mailing your Continuing Education Post Tests

Complete ALL hard copy unit post tests for the products you purchased in legible printing BEFORE your license expires. Mail is processed the same day it is received.

You may want to copy them BEFORE you mail them to the S.T.A.R.S. office to minimize mail delivery complications. They will NOT be returned to you unless you get a 75% or less. If you do NOT get a 75% or better after evaluation, the post tests will be sent back to you with the questions needing a new answer selection. After completing the questions, send them ALL back to the S.T.A.R.S. office for re-evaluation.

Be sure to use the CORRECT postage by having it weighed at the post office if it consists of more than 5 pages. Envelopes with INSUFFICIENT POSTAGE will be sent back to the participant and delay your post test evaluation and certificate creation.

I do NOT accept faxes since faxes fade over time and I need to keep my records for 3 years in case you would get audited by the Ohio Dept. of Health for some reason.

I do NOT accept scanned pages because I do NOT want you putting your private, personal information on the internet. I do NOT open attachments due to the threat of virus contamination that may jeopardize my web site and computerized databases.

Do NOT send your study media i.e. CD, DVD, booklets and/or books back to me.

United States Postal Service (USPS):

If you are using USPS for priority or express mailing, please keep your receipt with the tracking number in case of a problem with the delivery. Please mark the section for NO SIGNATURE REQUIRED for express mail and send it to my home address: Carolyn J. Frigmanski, MA, BSRT (R) 3134 Aldringham Road Toledo, Ohio 43606. The USPS does NOT deliver to my P.O. Box address. Please call to let me know I should be expecting it at 419 471-1973.

FedEx or United Parcel Service UPS:

If you are using these delivery services, please keep your receipts with the tracking number in case of a problem with the delivery. Please mark the section for NO SIGNATURE REQUIRED for express mail and send it to my home address: Carolyn J. Frigmanski, MA, BSRT(R) 3134 Aldringham Road Toledo, Ohio 43606. Please call to let me know I should be expecting it at 419 471-1973.

Thank you very much.



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo, Ohio 43606 419-471-1973

E-mail: info@xrayhomestudies.com

www.xrayhomestudies.com

How to renew your GXMO license in Ohio:

The Ohio Dept. of Health (ODH) requires a minimum of 12 continuing education credits (ceus) to be completed every 2 years (your biennium) BEFORE your license expires. You may do more than 12 ceus, but not less than 12 ceus, if you so choose. Ceus in excess of 12 cannot be carried over to the next biennium.

You will receive a hard copy renewal notice by mail from the ODH 60 days **BEFORE** your license expires. It is your responsibility to amend your personal information to the ODH whenever you change your name, address or place of employment as soon as possible by using the ODH website or contacting the ODH by telephone at 614 752-4319 for assistance. Failure to receive an ODH notice is not an acceptable reason for failing to renew on time. You can add completion of clinical modules to your GXMO license on the ODH web site.

Your ODH notice informs you that you may renew online or request a hard copy form from them. You must have your S.T.A.R.S. certificate(s) of completion immediately available when you renew since your course title(s), number of ceus, and ODH accreditation number(s) and date(s) of completion are printed on it.

You can renew immediately when you receive your notice or you have 30 days to complete the renewal process and payment to the ODH after your license expires. Online renewal requires your credit card for payment. If you chose hard copy renewal, you may submit a check or money order.

You and/or your employer can view and/or print your renewed license on line upon completion of the process. Problems that exist with renewal should be addressed to the ODH by calling for assistance.

S.T.A.R.S. personnel are **NOT** responsible for your renewal. Please direct any questions or needed assistance with renewal to the ODH personnel.

GXMOs must notify the ODH in writing within 30 days of any changes in the physician providing direct supervision. If your scope of practice changes (e.g. chiropractic to podiatric), a competency form must be completed and submitted to the ODH.

You may also want to check the ODH web site periodically for changes that may have occurred during your biennium and to share this information with your co-workers and/or administrative staff members.

The ODH website is: http://www.odh.ohio.gov/odhPrograms/rp/rlic/ristatus.aspx

Email is: BRadiation@odh.ohio.gov

Thank you very much.



Specialized Topics in Areas of Radiologic Sciences

P.O. Box 2931 Toledo, Ohio 43606

Phone: 419-471-1973

Website: www.xrayhomestudies.com

T	Jni	t 3	6

Title: Radiographic Imaging & Exposure by Terri L. Fauber

Please complete the answer sheet at the conclusion of this post test and return to the S.T.A.R.S. office.

Chapter 1

1.	The crystals Dr.	Wilhelm	Conrad	Roentgen	used to	produce	the	first x-ray	v were
----	------------------	---------	--------	----------	---------	---------	-----	-------------	--------

a. calcium tungstate

c. lanthanum oxybromide

b. barium platinocyanide

d. gadolinium oxysulfide

2. The instantaneous production of light resulting from some type of energy is called

a. luminescence

c. fluorescence

b. phosphorescence

d. scintillation

3. The world's first radiograph was a hand image of

a. Dr. Wilhelm Conrad Roentgen

c. Marie Curie

b. Thomas Edison

d. Anna Bertha Roentgen

4. The state that attempted to enact legislation to ban x-ray producing devices in opera glasses was

a. New Jersey

c. New York

b. Washington, DC

d. New Mexico

5. The first indications of the damaging biologic effects from x-rays were

a. tumors & sterility

c. mental challenges & fatigue

b. redness & burns of the skin

d. soft tissue tumors and anemia

6. The most energetic electromagnetic radiation on the electromagnetic spectrum is

a. cosmic radiation

c. ultraviolet

b. infrared radiation

d. gamma rays

7. Constant velocity is

a. the speed of sound

c. 186,000 miles per second

b. the frequency of electricity

d. the cycles of all electromagnetic energy

8. The maximum energy of a photon is expressed as

a. the speed of sound

c. the speed of light

b. wavelength

d. kilovoltage peak (kVp)

9.	X-ray:	s can cause chemical and bio	ological changes	in living tissue because of the process of
20	. a.	ionization	c.	scatter
	b.	reduction	d.	oxidation
Cha	apter i	2		
10.	The co	oiled tungsten wire that prov	ides electrons du	ring x-ray production is called the
	a.	cathode	C.	anode
	b.	filament	d.	focusing cup
11.	The re	volutions per minute of rota	ting anodes is	
	a.	1,000 to 2,000	c.	3,000 to 10,000
	· b.	500 to 2,500	d.	300 to 1,000
12.	Tungs	ten is used in the target beca	use it possesses	a
	a.	high atomic number		lovy boot store so populity
	а. b.	low atomic number		low heat storage capacity limited capacity to produce scatter
.13.	The p	hysical area of the target tha		y electrons during x-ray production is
	1	•	~ 10 ~ 0 111 0 11 0 10 0 0 0 0	J elseword daring h ray production is
	a. h	focal track focusing cup		stator focal spot
14.			•	interact with K-shell electrons is called
	a.	characteristic	c.	bremsstrahlung
	b .	x-ray emission	d.	Compton effect
15.	The po	ortion of the x-ray circuit tha	t operates on 3 to	5 amps and about 10 volts is called the
	a.	high voltage circuit	c.	low voltage circuit
	b.	filament circuit	d.	rectification system
16. I	Electro	on flow from cathode to anot	de in the x-ray tu	be is called
. : •	a.	actual focal spot	C.	effective focal spot
	b.	high voltage flow	d.	tube current
17. I	Beam j	penetrability increases as		
	a.	mA increases	. C.	kVp increases
	b.	filter thickness increases	d.	kVp decreases
18. 7	Гһе ех	posure factor that produces	or controls the m	umber of x-rays is
•	a.	distance	C.	kVp
		mA	•	distance

19. The ar	node angle typically ranges from		
. a.	10 to 15 degrees	C	5 to 15 degrees
	6 to 18 degrees		6 to 20 degrees
		-	
20. Thick	er anatomic body parts should be position	ed	at the
a.	cathode end of the x-ray tube	C.	anode end of the x-ray tube
ь.	head of the x-ray table		transverse axis of a grid
21. Total	filtration for x-ray tubes operating above	70	kVp in the United States is
a.	1.5 mm of aluminum	C.	2.5 mm of aluminum or its equivalent
	2 mm of aluminum or its equivalent		3.5 mm of aluminum
	e of the rotor to rotate the anode will caus	e	
	pitting		uneven x-ray beam distribution
b.	melting of the focal tract	d.	damage to the cathode structure
Chapter 3 23. The re		am	as it passes through anatomic tissue is called
a.	attenuation	c.	absorption
b.	differential absorption	d.	braking radiation
24. The pr	ocess of creating photoelectrons occurs w	/ith	
a.	Compton effect	C.	exit radiation
b.	scatter	d.	photoelectric effect
25. The co	mbination of transmitted and scattered ra	dia	tion is called
			primary radiation
b.	absorbed radiation	d.	exit radiation
26. The im	naging equipment involving the interaction	n o	f radiation with an input phosphor is called
			digital radiography
b.	image intensification	d.	computed radiography
27. A spec	ialized image receptor that can produce a	CO	mputerized radiographic image is utilized in
a.	10		xerography
b.	image intensification	d.	digital radiography
Chapter 4 28. The vis	sibility of recorded detail is also referred t	to a	S

a. geometric propertiesb. sensitometric properties

c. photographic propertiesd. densitometric properties

29. The c	controlling factors for radiographic density	y are
	mA & time kVp & filtration	c. time & distance d. kVp & distance
30. If rep	eat radiographs are necessary because of a	a minimal density error, the mAs should be adjusted
a. b.	by a factor of 3 by a factor of 2	c. by a factor of 1 and ½ d. by a factor of 4
31. The p	percentage of kVp that must be used to ma	intain density when mAs is doubled is
a. b.	20 10	c. 30 d. 15
32. The f	ormula to use to maintain consistent dens	ity when the source-to-image distance is changed is
	density maintenance formula reciprocity law	c. inverse square law formula d. electrostatic laws
33. When	switching from non-grid to a 12:1 ratio g	rid, mAs needs to be increased by a factor of
a.	2 b. 4	c. 5 d. 3
34. Small	er field size/increased collimation will af	fect density in what way?
	increase radiographic density decrease radiographic density	c. increase scatter d. decrease contrast
35. mAs s	should be doubled when part thickness inc	creases by
a.	2 cm b. 10 cm	c. 4 cm d. 6 cm
36. To pro	oduce a more uniform density on varying	body thicknesses, the radiographer can use
a. b.	compensatory filters higher kVp	c. the reciprocity law d. cathode heel effect
37. A radi	ograph possessing few densities with great	at differences among them is described as
a. b.	long- scale contrast average contrast	c. low contrast d. short-scale contrast
38. The co	ontrolling factor of contrast is	
a.	mA b. time	c. kVp d. distance
39. Short-	scale contrast can be created by utilizing	
a. b.	low kVp a low grid ratio	c. high mA d. a larger collimator size

40. The	production of scatter radiation	increases when	ever	
	. kVp decreases . anatomic tissue increases		part thickness do	
41. The	type of contrast agent that produced	duces less radio	graphic density th	an the adjacent tissue is
a	. negative b. bariu	m sulfate c	air	d. positive
42. The	mage receptor most sensitive	to scatter radiat	ion is used in	
a b	a film-screen system image intensification		digital imaging fluoroscopy	
43. Adul	t mAs exposure factors can be	adjusted for pe	diatric exams und	ler age 5 by utilizing
	25% of the adult technique 10% of the adult technique		50% of the adult 35% of the adult	
44. Dry _I	laster casts require an increaso	e in mAs of		
a	3x b. 2x	C.	4 x	d. 1 1/2x
45. The c	lisease considered to be a dest	ructive conditio	n is	
	pleural effusion cirrhosis		congestive heart emphysema	failure
46. Soft 1	issue technique would be nece	essary to use for	patients experien	cing
a. b.	croup tuberculosis		subcutaneous em hemorrhage	physema
Chapter 47. The d	5 istinctness or sharpness of stru	ıctural lines tha	t make up the radi	iographic image is called
	geometric properties photographic properties	*	recorded detail resolution	
48. The a	verage focal spot dimensions a	are in which of	the following rang	ges?
a.	1 to 2 cm b. 0.5 to	1.2 mm c.	2 to 3 mm	d. 1 to 3 cm
49. Increa	sing SID and decreasing OID	will result in _		in an image.
	increasing the amount of sha less structural detail	_	reducing patient of decreasing the an	doses nount of unsharpness
50. Increa	sing the speed of an intensifyi	ng film-screen :	system will	
a. b.	decrease recorded detail increase patient doses		increase recorded improve quality c	

51. The v	isual acuity of the human eye can discern	ap	proximately
	2 line pairs per millimeter 3 line pairs per millimeter		. 5 line pairs per millimeter . 4 line pairs per millimeter
52. The n	nost detrimental factor to maximizing reco	ord	ed detail is
	incorrect processing incorrect distance selection		improper grid selection patient motion
53. To pro	oduce a magnified image, the radiographe	er n	nust
	decrease SID increase OID		use screens use the anode heel effect
54. When	utilizing a 72" SID and a 60" SOD, the n	nag	mification factor is
a.	2 b. 1.2	C.	d. 2.5
55. To mi	nimize shape distortion, the radiographer	mı	st keep the tube, film and patient
	perpendicular to each other horizontal to the table		as close as possible properly aligned
56. Increa	sing the angle on the central ray will crea	te a	
	decrease in density increase detail		increase in density decrease distortion
Chapter 6 57. Increa	sed scatter will be produced whenever		
	patient thickness is decreased collimation is small		tissue volume increases a destructive disease is present in a patient
58. As bea	m restriction or collimation increases, the	e pa	atient dose will
	increase decrease by at least 50%		will not be affected decrease
59. Signif	icant collimation requires an increase in r	mA	s of
	20 to 40% 30 to 50%		10 to 15% 60% or more
60. The be	am-restricting device that utilizes a pair o	of a	djustable lead shutters is a
a. b.	collimator cylinder cone		diaphragm flared extension cone

61. Another term for an automatic collimator is	a/an device.
a. adjustableb. electronic	c. positive beam-limiting d. self-regulating
62. A grid should be utilized whenever an anator	mic body part measures
a. 15 cm or moreb. 10 cm or more	c. 5 to 9 cm d. 8 to 9 cm
63. The grid type that has a convergence point is	the control of the co
a. parallel b. cross-hatched	c. wafer d. linear focused
64. The grid conversion factor when 100 mAs is	used with a grid and 5 mAs is used non-grid is
a. 10 b. 15	c. 500 d. 20
65. Misalignment of the grid in relationship to the	ne primary x-ray beam results in
a. grid cutoffb. increased contrast improvement	c. increased density d. reduced patient dose
66. The grid ratios typically used in diagnostic ra	diography are
a. 6:1 & 8:1 b. 5:1 & 10:1	c. 8:1 & 12:1 d. 8:1 & 16:1
67. Digital imaging systems, when compared to	film-screen imaging systems, are more responsive to
a. high energy photonsb. very low energy photons	c. scattered photons d. photoelectrons
Chapter 7 68. The most common film-screen image receptor	or utilized today is
a. digital b. single-emulsion film	c. low-dose film d. double-emulsion film
69. Flat silver halide crystals are called	
a. T-grain b. crystalline matrix	c. horizontal d. conventional
70. The percentage of silver bromide in film em	ulsions used by manufacturers today is
a. 75 to 90 b. 50 to 75	c. 90 to 99 d. 60 to 80
71. The layer of film construction that actually cr	eates the image is called the

c. base

d. photographic layer

a. supercoatb. emulsion

72. The i	ocal point for the development of the la	tent image centers is fou	ınd in the
a b		c. silver bromide crd. sensitivity specks	
73. The t	erm used to describe an image being rec	orded on a film by refle	cted light is
a. b.	halation image intensification	c. anti-crossover d. photo-reception	
74. A filr	n-screen system that is not spectrally ma	atched will produce dens	sities that are considered
a. b.	optimal diagnostic	c. suboptimal d. extremely desirab	ble
75. The e	ffect of crossover creates images with		
	maximized detail a loss of recorded detail	c. less distortion d. better contrast sca	hles
76. The e	mission of light from the screens when	stimulated by x-rays is c	alled
a. b.	T	c. phosphorescence d. luminescence	
77. Lanth	anum oxybromide is a crystal in the cate	egory of	·
a. b.	blue emitting rare earth conventional blue emitting	c. green emitting rar d. ultraviolet emittin	
78. The in	ntensification factor when using 5 mAs v	with screens and 40 mAs	s without screens is
. a.	10 b. 20	c. 8	d. 0.8
79. A scre	een is considered fast when its relative s	peed is at least	
a.	200 or more b. 400 or more	c. 150 to 300	d. 80 to 200
80. A ligh	t-absorbing dye added to the phosphor l	ayer will	
a. b.	create fog on the resultant image improve the scale of contrast	c. increase film speed. improve recorded	
81. A syne	onym for image noise is		
a. b.	quantum mottle photon interaction	c. screen crossover d. speed conversion	
82. The qu	uality control device used to check film-	screen contact is called a	a

c. wire mesh test toold. ionization chamber

a. spinning topb. penetrometer

	83. The fu	nction of lead foil in the back of the ca	ssette is to	
	. a. b.	prevent crossover absorb back scatter	c. prevent phosphorescence d. absorb light photons	
	84. A pho	tostimulable phosphor is used in		
	A	digital radiography	c. computed tomography	
	a. b.	xerography	d. computed radiography (CR)	
	85. Bariur	n fluorohalide crystals coated with euro	opium are found in	
	a.	rare earth screens	c. imaging plates	
•	b .	input phosphors	d. super fast screen speeds	
· . : .	86. To pre	vent fading, computed radiography (Cl	R) plates should be processed within	
	·. a.	1 hour of exposure	c. 2 hours of exposure	i
	ь.	30 minutes of exposure	d. 3 hours of exposure	
	87. X-rays	are absorbed by a layer of cesium iodi	de in	
	a. b.	a direct conversion system an indirect conversion system	c. all computed radiography systems d. all fluoroscopy systems	
	Chapter 8 88. The te	B rm used to describe an invisible image	on exposed film is called	
	a.	manifest	c. crossover	
	b.	photosensitivity	d. latent	
	89. The ch	nemical that is fast -reducing and produ	ces gray densities on film is	
:	a.	hydroquinone	c. phenidone	
	b.	sodium sulfite	d. potassium bromide	
	90. The ch	nemical added to the processing solution	ns to accommodate the roller transport system is	
	a.	glutaraldehyde	c. hydroquinone	
	b.	sodium sulfite	d. ammonium thiosulfate	
	91. The ar	nmonia odor emanating from processin	ng solutions is attributed to the	
	a.	preservative	c. acidifier	
	b.	fixing agent	d. hardener	
. !	92. Stainir	ng or fading of the permanent or archiva	al image results from	
	a.	improper washing	c. contamination of the solutions	
	b .	improper development	d. high humidity in the file room	

93. The ch	emical that actua	lly stops developmer	nt is			
a. b.	sodium sulfite chrome alum	7. 		potassium bromi acetic acid	de	
94. To enh	ance archival qua	lity, finished radiogr	raphs :	should maintain		
·	20 to 25% water 0% water			10 to 15% water 5 to 8% water		
95. The rol	ller assembly that	moves film from on	e tank	to another is calle	d the	
	turnaround crossover			deep or planetary entrance		
96. The rec	luction of chemic	al strength as a resul	t of ex	posure to air is cal	lled	
b.	aerial oxidation exhaustion			aerial reduction evaporation		
97. Repleni	ishment pumps ar	e activated by the			2	
	feed tray re-circulation pur	nps		entrance rollers micro-switch		
98. Temper	rature control in n	ewer automatic proc	essors	is provided by the	;	
b.	mixing valves external thermost			immersion heater dryer assembly		
99. A brown	n or thiosulfate st	ain results from				
	inadequate washi exhausted develo			contamination of sinadequate fixatio		
100. One of	the causes of dic	hroic fog is				· ·
and the second s	mproper washing ınder-replenishm			over-replenishmer insufficient harder		
101. Unused	d film should be s	tored in a relative hu	ımidit	y range of		
a. 2	20 to 30%	b. 40 to 60%	c.	50 to 70%	d. 0 to 15%	
102. Archiv	al storage temper	atures in Fahrenheit	degree	es for radiographic	film ranges be	tween
•	50 to 60	b. 60 to 70		60 to 80	d. 55 to 75	•
103. A com	monly used filter	in darkrooms storing	g blue-	sensitive film is ca	ılled	
	Kodak Wratten 6E			Kodak GBX		

d. Kodak GB

b. Kodak Wratten 15 watt

104. Quali	ty control for autor	natic processors requi	res	replenishment rate	che	cks every	
· a.	day	b. week	C.	quarter	d.	year	
105. Quali	ty control for autor	natic processors requi	res	developer and fixer	рН	checks every	
*********** a.	day	b. week	C.	quarter	d.	year	
106. The m	nost efficient metho	od of silver recovery is	3				
	metallic replacements	ent		silver filtration electrolytic			
107. Bendi	ng or kinking of fi	lm will create an artifa	ict i	identifiable as			
	an abrasion a half-moon mark			an area of fog a minus density			
108. Foreig	n objects are visua	lized on radiographic	filn	n as a			
a. b.	positive-density ar radiolucent artifac	tifact		minus-density artif an area of fog	act		
Chapter 9 109. A devi	ice used to produce	e consistent step-wedg	e de	ensities by eliminat	ing	the x-ray unit is	a
	sensitometer penetrometer			densitometer wire mesh			
110. For dia	agnostic ability, ph	ysicians need to read o	opti	cal densities in the	rang	ges of	•
a.	1.0 to 2.5	b. 1.5 to 3.5	c.	1.5 to 3	d. (0.25 to 2.5	
111. A doul	bling or halving ch	ange occurs in density	wh	nen the log relative	ехр	osure changes b	У
a. 1	1.0	b. 0.3	C.	1.5	d. :	2.0 or more	* * * * * * * * * * * * * * * * * * *
112. The reg	gion in which a ch	ange in exposure inten	sity	no longer affects o	ptic	cal density is cal	lled
	he shoulder O _{max}		c. :	straight-line toe			
113. When o	comparing two film	n products, the lower s	spec	ed exposure point w	'ill i	ndicate the film	with
	slower speed petter resolution			faster speed improved contrast s	cale	es	
114. When o	comparing two film	n products, the steeper	the	slope will indicate	the	film with	
	nigher contrast setter resolution			slower speed less vulnerability to	fog	ging	

115. Naı	row latitude film products possess	
a b		c. higher densitiesd. higher contrast scales
116. To	achieve maximum contrast, optical densit	ies on a sensitometric curve must lie within
a		c. the toe region d. the speed point region
117. The	term used to describe the shades of gray	in a digital image is called the
a b	. average gradient . shoulder range	c. dynamic range d. contrast enhancement scale
Chapter 118. Bas	· 10 eline kVp should have 2 kVp added for ev	very increase of patient tissue of
a	. 1 cm b. 3 cm	c. 2 cm d. 4 cm
119. Vai	iable kVp/Fixed mAs technique charts are	recommended for use with
	. abdomen studies . chest radiography	c. small extremities d. gastrointestinal tract examinations
120. The	kVp value that is high enough to ensure p	penetration of the part is called
a b		c. diagnostic d. useful
121. For	exposure technique charts to be effective,	radiographic equipment must be working
a b		c. within 5-10% of the calibration parameters d. within normal limits
122. Dev	velopment and utilization of exposure tech	nique charts serve to
	. provide consistent quality radiographs . accommodate all pathologies	c. accommodate all patient sizes d. provide no personnel errors
Chapter 123. The	11 electronic device that converts visible lig	ht energy into electrical energy is
.a	. an ionization chamber . AEC detector	c. photomultiplier tube (PM) d. input phosphor
124. Mo	st of today's automatic exposure systems ((AEC) utilize

c. one AEC chamber

d. excessive kVp to penetrate all parts

a. 10 photomultiplier tubesb. ionization chambers

	125. A cr	itical factor for radiographers to reme	mber when using an AEC sy	stem is to
d .			- •	•
	a.	L - L 3	 c. disregard using gon 	adal shielding
	ъ.	disregard collimation	d. select all chambers	for all exams
				- WIL WARREN
	126. Ano	ther critical factor for radiographers to	remember when using an A	EC system is to
	a.	adjust time for disease processes	c, disregard SID	
• . •		always use high kVp	d. select the correct ch	amhar
			poreof the collect Cli	
	127. Dens	sity selections with AEC systems can	be adjusted in increments of	•
	a.	10 to 20%	c. 30 to 50%	
	. b.	25 to 50%	d. 10 kVp	
		•		
	128. The	device that prevents the x-ray tube fro	m reaching or exceeding its	heat loading capacity is
	a.	the back up timer	c. the rheostat	
		the ionization chamber	d. the rectifiers	
	•		en manifeld	•
	129. Radio	ographers should note the mAs readou	ut display so that they can	
	a.	perform all quality control tests	c. make adjustments fo	manual tacketers
	b.		d. tell the patient the d	or manual techniques
•	3.	aromina without accommended	d. ten me panem me u	ose.
	130. A sys	stem using a preprogrammed set of ex	posure factors with visual ar	natomic buttons is called
•	a.	anatomic AEC system	c. preprogrammed AE	ariatowa .
	b.	anatomic adjustment system	d. anatomic programm	9
		and the state of t	d. anatomic programm	ıng
	131. The s	ensitivity of CR systems is approxima	ately equivalent to a film-scr	een speed system of
•		300 b. 200	400	
	a.	300 b. 200	c. 400	d. 100
· · . · · ·				
	Chapter 1			
	132. Digita	al images are constructed from		
in the second	n de la companya de La companya de la co			
		numerical data	c. megabytes	
	b .	electronic charges	d. x-rays	
	•			
• • •	133. The p	rocess of manipulating the image qua	lity such as contrast, subtrac	tion, etc. is called
	a.	pixel processing	c. postprocessing image	anhonoomost
•		matrix manipulation		eminicement
	υ.	many manymanon	d. pixel preprocessing	
	134. The c	reation of metallic silver in radiograpl	nic images would correlate to	
	a.	pixel processing	a machana a main a s	1.
			c. postprocessing image	ennancement
,	υ.	matrix manipulation	d. pixel formation	

135. Incre	easing the matrix size and pixel number will	
	improve image quality	c. distract from image quality
ь.	minimize image manipulation ability	d. decrease image quality
136. The	analog image is converted into a digital image f	or computer processing by the
a.	charged-coupled device	c. imaging plate
b.	reader unit	d. processing unit
137. One	of the disadvantages of computed radiography ((CR) in comparison to film-screen systems is
a.	the potential for increased patient exposure	c. inability to see fine recorded detail
ь.	lack of image sizes	d. inability to use with contrast agents
138. The	digital composite of the varying x-ray intensities	s existing in the patient is identified by the
a.	analog-to-digital converter	c. scintillating device
b.	charged-coupled device	d. image matrix
139. The g	greatest detective quantum efficiency (DQE) is	found in
a.	film-screen systems	c. flat panel direct capture systems
. b.	computed tomography	d. automatic exposure systems
140. A po	stprocessing technique that improves the visibil	ity of small, high contract structures is called
a.	subtraction	c. contrast enhancement
b.	edge enhancement	d. black/white reversal
%		



Specialized Topics in Areas of Radiologic Sciences P.O. Box 2931 Toledo, Ohio 43606

Phone: 419-471-1973

Website: www.xrayhomestudies.com

from S.T.A.R.S. wh	en you receive a	nation so that you can obtain 75% or higher score.	n a signed certificate from an official (Please print)	
Name	,			
Address				***************************************
City		State	Zip Code	
Social Security Number		Date	v.	

Answer Sheet for: Unit 36 Radiographic Imaging & Exposure by Terri L. Fauber

Please place your lettered selection for each question in the respective box and return ONLY this post test sheet to S.T.A.R.S.

1.	11.	21.	31.	41.	51.	61.
2.	12.	22.	32.	42.	52.	62.
3.	13.	23.	33.	43.′	53.	63.
4.	14.	24.	34.	44.	54.	64.
5.	15.	25.	35.	45.	55.	65.
6.	16.	26.	36.	46.	56.	66.
7.	17.	27.	37.	47.	57.	67.
8.	18.	28.	38.	48.	58.	68.
9.	19.	29.	39.	49.	59.	69.
10.	20.	30.	40.	50.	60.	70.

Please continue on the reverse side.

Answer Sheet for: Unit 36 Radiographic Imaging & Exposure by Terri L. Fauber

Please place your lettered selection for each question in the respective box and return ONLY this post test sheet to S.T.A.R.S.

71.	81.	91.	101.	111.	121.	131.
72.	82.	92.	102.	112.	122.	132.
73.	83.	93.	103.	113.	123.	133.
74.	84.	94.	104.	114.	124.	134.
75.	85.	95.	105.	115.	125.	135.
76.	86.	96.	106.	116.	126.	136.
77.	87.	97.	107.	117.	127.	137.
78.	88.	98.	108.	118.	128.	138.
79.	89.	99.	109.	119.	129.	139.
80.	90.	100.	110.	120.	130.	140.